

PERFORMANCE CHARACTERISTICS

A) Internal Evaluation:

- Accuracy: In an internal study **Electra™ AFP** was evaluated against commercially available licensed kit with 90 random clinical samples, & **Electra™ AFP** has demonstrated 100% clinical correlation with the commercially available licensed kit.
- Precision: **Electra™ AFP** was evaluated with licensed external Quality controls for Precision Studies & following is the data:

Controls	No. of testings	Mean Control values with Electra™ AFP	Coefficient of Variation (CV)
Level 1	10	8.32	4.23
Level 2	10	69.19	2.76
Level 3	10	196.02	1.44

B) External Evaluation:

Electra™ AFP CLIA has been evaluated by a NABL accredited lab against their reference method. In this evaluation **Electra™ AFP** has demonstrated 100% correlation with the reference method.

*Data file: Zephyr Biomedicals (A Division of Tulip Diagnostics Pvt. Ltd).

Important Note:












- The **Electra™ AFP** assay is a temperature sensitive assay. The best temperature condition for this assay is from 18°C to 25°C.
- The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.
- It is recommended to use the multiple channel pipettes to avoid time effect. A full plate of 96 wells may be used if automated pipette is available.
- Duplication of Standards & samples is not mandatory but may provide information on reproducibility & application errors.

LIMITATIONS OF THE ASSAY

(1). As with all diagnostic tests, a definite clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated. (2). The activity of the enzyme used is temperature-dependent and the RLU values may vary. The higher the room temperature (+18°C to +25°C) during substrate incubation, the greater will be the RLU values. Corresponding variations apply also to the incubation times. However, the Standards are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result. (3). Adaptation of this assay for use with automated sample processors and other liquid handling devices, in whole or in part, may yield differences in test results from those obtained using the manual procedure. It is the responsibility of each laboratory to validate that their automated procedure yields test results within acceptable limits. (4). Insufficient washing (e.g., less than 5 wash cycles, too small wash buffer volumes, or shortened reaction times) can lead to incorrect results.

BIBLIOGRAPHY

(1). Rongen HA, Hoetelmans RM, Bult A, van Bennekom WP. Chemiluminescence and immunoassays. **J Pharm Biomed Anal** 1994 Apr;12(4):433-62. (2). Koszegi T. Immunoluminometric detection of human prolactin. **J Biochem Biophys Methods** 2002 Oct;53(13):157-64. (3). Roda A, Simoni P, Mirasoli M, Baraldini M, Violante FS. Development of a chemiluminescent enzyme immunoassay for urinary 1-hydroxypyrene. **Anal Bioanal Chem** 2002 Apr;372(7-8):1481-1486. (4). Laffin RJ, et al. Hybritech total and free prostate-specific antigen assays developed for the Beckman Coulter access automated chemiluminescent immunoassay system: a multicenter evaluation of analytical performance. (5). **Clin Chem** 2001 Jan;47(1):129-32. (6). Koszegi T, Walker W.H.C. Introduction: An Approach to Immunoassay. **Clin. Chem.** 1977; 23: 384. (7). Kirkegaard C., Friis T. and Siersback-Nielsen K. **Acta Endocrinol.** 1974; 77: 71. (8). Wisdom G.B. Enzyme-Immunoassay. **Clin. Chem.** 1976; 22: 1243. (9). Hoffenberg R. **Medicine** 1978; 8: 392. (10). Lieblisch J., Utiger R.D. **J. Clin. Invest.** 1972; 51: 1939.

 Temperature Limitation	 Consult Instructions for use	 Date of Manufacture	 Batch Number / Lot Number
 Manufacturer	 In vitro Diagnostic Medical Device	 This side up	 Contains sufficient for <n> tests
 Use by	 Catalogue Number	 Do not reuse	

Manufactured by:

Zephyr Biomedicals

A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa- 403 722, INDIA.

Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.

An ISO 13485
Certified Company

0422/VER-01

electra™
●●●●●●●● AFP

Chemiluminescence Assay for the Quantitative Determination of Alpha-fetoprotein (AFP) in Human Serum.
FOR IN VITRO DIAGNOSTIC USE ONLY
Store at 2°C to 8°C

INTENDED USE

ELECTRA™ AFP CLIA test is intended for the quantitative determination of Alpha-fetoprotein (AFP) in human serum. For In Vitro Diagnostic Use only.

INTRODUCTION

Alpha-fetoprotein (AFP) is normally produced during fetal and neonatal development by the liver, yolk sac, and in small concentrations by the gastrointestinal tract. After birth, serum AFP concentrations decrease rapidly, and by the second year of life and thereafter only trace amounts are normally detected in serum. Elevation of serum AFP to abnormally high values occurs in several malignant diseases, most notably nonseminomatous testicular cancer and primary hepatocellular carcinoma. In the case of nonseminomatous testicular cancer, a direct relationship has been observed between the incidence of elevated AFP levels and the stage of disease. Elevated AFP levels have also been observed in patients diagnosed with seminoma with nonseminomatous elements, but not in patients with pure seminoma. In addition, elevated serum AFP concentrations have been measured in patients with other noncancerous diseases, including ataxia telangiectasia, hereditary tyrosinemia, neonatal hyperbilirubinemia, acute viral hepatitis, chronic active hepatitis, and cirrhosis. Elevated serum AFP concentrations are also observed in pregnant women. Therefore, AFP measurements are not recommended for use as a screening procedure to detect the presence of cancer in the general population.

PRINCIPLE

ELECTRA™ AFP Quantitative CLIA assay is for use on **ELECTRA™** analyzer. **ELECTRA™ AFP CLIA** works on the principle of chemiluminescence wherein light is produced by a chemical reaction from a substance as it returns from an electronically excited state to the ground state. When catalysed by HRP, the oxidation of luminol by hydrogen peroxide produces an electronically excited form of 3-aminophthalate which on relaxation emits light with maximum intensity at $\lambda=425\text{nm}$. In **ELECTRA™ AFP**, the assay system utilizes one Anti-AFP antibody for solid phase (microtiter wells) immobilization and another mouse monoclonal Anti-AFP antibody in the antibody-enzyme (horseradish peroxidase) conjugate solution. The test specimen (serum) is added to the AFP antibody coated microtiter wells and incubated with the zero buffer. If human AFP is present in the specimen, it will combine with the antibody in the well. The well is then washed to remove any residual test specimen, and AFP antibody labeled with horseradish peroxidase (conjugate) is added. The conjugate will bind immunologically to the AFP on the wells, resulting in the AFP molecules being sandwiched between the solid phase and enzyme-linked antibodies. After incubation at room temperature, the wells are washed with wash buffer and bound enzyme is detected by adding the chemiluminescent substrate. The bound enzyme converts substrate to a reaction product that emits a photon of light. Chemiluminescence is measured in Relative Light Units (RLU). The amount of light emitted is proportional to the amount of enzyme present and is directly related to the amount of AFP in the sample. By reference to a series of AFP standards assayed in the same way, the concentration of AFP in the unknown sample is quantified.

MATERIALS & COMPONENTS

Materials provided with the test kits:

- Coated Microwells: Microwells coated with Anti-AFP antibody.
- AFP Zero Buffer. Ready to use.
- AFP Enzyme Conjugate. Ready to use.
- Substrate A: Chemiluminescent substrate containing enhanced luminol solution.
- Substrate B: Chemiluminescent substrate containing stabilized peroxide solution.
- AFP Standard set of 6 standards labeled as A to F in liquid form. Ready to use. For standard Concentrations refer vial label.
- Wash Buffer Concentrate (20X).

Materials required but not provided:

- Precision pipettes: 10-100 μl , 20-200 μl , 100-1000 μl
- Disposable pipette tips
- Distilled water
- Disposable Gloves
- ELECTRA™** Analyzer

STORAGE AND STABILITY

- ELECTRA™ AFP** kit is stable at 2-8°C up to the expiry date printed on the label.
- Coated micro-wells should be used within one month of opening the pouch. Once opened, the pouch must be sealed.

properly to protect from moisture. In case the desiccant pouch changes color from blue to white, the strips should not be used.

3. Diluted wash buffer is stable up to one week at 2-8°C.
4. Working Substrate (A+B) must be used immediately.

SPECIMEN COLLECTION

1. Collect blood specimen by venipuncture according to the standard procedure.
2. Only serum should be used.
3. Avoid grossly hemolytic, lipemic or turbid samples.
4. Preferably use fresh samples. However, specimens can be stored up to 48 hours at 2-8°C, for short duration.
5. For longer storage, specimens can be frozen at -20°C. Thawed samples must be mixed prior to testing.
6. Do not heat inactivate before use.
7. Specimen containing precipitate or particulate matter should be clarified by centrifugation prior to use.
8. Specimen should be free from particulate matter and microbial contamination

PRECAUTIONS

1. Bring all reagents and specimen to room temperature before use.
2. Do not pipette any material by mouth.
3. Do not eat, drink or smoke in the area where testing is done.
4. Use protective clothing and wear gloves when handling samples.
5. Use absorbent sheet to cover the working area.
6. Immediately clean up any spills with sodium hypochlorite.
7. All specimens and standards should be considered potentially infectious and discarded appropriately.
8. Neutralize acid containing waste before adding hypochlorite.
9. Do not use kit after the expiry date.
10. Do not mix components of one kit with another.
11. Always use new tip for each specimen and reagent.
12. Do not allow liquid from one well to mix with other wells.
13. Do not let the strips dry in between the steps.


REAGENT PREPARATION

- All reagents should be brought to room temperature (18-25°C) and mixed by gently inverting or swirling prior to use. Do not induce foaming.
- Dilute wash buffer 20 times (for example add 5ml concentrated buffer to 95 ml distilled or deionized water). Mix well before use.
- Prepare a Working Substrate by Mixing Substrate A and Substrate B in equal volume (1:1 ratio) before addition to the micro-wells

No. of Strips	1	2	3	4	5	6	7	8	9	10	11	12
Substrate- A µl	250	450	650	850	1050	1250	1450	1650	1850	2050	2250	2450
Substrate- B µl	250	450	650	850	1050	1250	1450	1650	1850	2050	2250	2450

TEST PROCEDURE

1. Secure the desired number of coated wells in the holder. Dispense **25µl** of Standards & serums into the appropriate wells.
2. Dispense **100 µl** of zero buffer into each well. Incubate at room temperature (18-25° C) for **15 minutes**.
3. After incubation, empty the microtitre wells and wash the plate 5 times with 350µl of diluted wash buffer. Strike the microtitre plate sharply onto absorbent paper towel to remove all residual droplets.
4. Dispense **100µl** of Enzyme Conjugate into each well. Incubate at room temperature (18-25°C) for **15 minutes**.
5. After incubation, empty the microtitre wells and wash the plate 5 times with 350µl of diluted wash buffer. Strike the microtitre plate sharply onto absorbent paper towel to remove all residual droplets.
6. Add **50 µl** of working Substrate (A+B) in all the micro-wells. Keep away from direct light while adding the substrate.
7. Cover the **ELECTRA™** microplate and incubate for **10 minutes** at room temperature (**18-25°C**) in dark.
8. Read the **ELECTRA™** micro-plate exactly at **10 minutes** in **ELECTRA™ Analyzer**.



CALCULATIONS


Construct a standard curve by plotting the RLU obtained from each reference standards against its concentration in ng/ml on the graph paper, with RLU values on the vertical or Y axis and concentrations on the horizontal or X axis. Use the RLU values for each specimen to determine the corresponding concentration of AFP in ng/ml from the standard curve. Any diluted specimens must be corrected by the appropriate dilution factor.

Example of Standard curve

Results of a typical standard run with RLU's shown in the Y axis against AFP concentrations shown in the X axis.

Suggest: Use 4-Parameter Standard curve to calculate sample values.

AFP Values (ng/ml)	RLU's
A	462
B	94120
C	502859
D	1448474
E	4214755
F	6655116



This Standard curve is for the purpose of illustration only and should not be used to calculate unknowns. Each user should obtain their own Standard curve and data.

Expected Ranges of values

In high-risk patients, AFP values between 100 and 350 ng/ml suggest a diagnosis of hepatocellular carcinoma, and levels over 350 ng/ml usually indicate the disease. Approximately 97% of the healthy subjects have AFP levels less than 8.5 ng/ml. It is recommended that each laboratory establish its own normal range. The minimum detectable concentration of AFP by this assay is estimated to be 2.0 ng/ml.