

PERFORMANCE CHARACTERISTICS

A) Internal Evaluation:

- Accuracy: In an internal study **ELECTRA™ Ferritin kit** was evaluated against commercially available licensed kit with 90 random clinical samples & **ELECTRA™ Ferritin** has demonstrated 95% clinical correlation with the commercially available licensed kit.
- Precision: **ELECTRA™ Ferritin kit** was evaluated with licensed external Quality controls for Precision Studies & following is the data:

Controls	No. of testings	Mean Control values with ELECTRA™ Ferritin	Coefficient of Variation (CV)
Level 1	10	75.44	3.36
Level 2	10	146.53	3.84
Level 3	10	324.40	5.21

B) External Evaluation:

ELECTRA™ Ferritin CLIA has been evaluated by a NABL accredited lab against their reference method. In this evaluation **ELECTRA™ Ferritin CLIA** has demonstrated 95% correlation with the reference method.

*Data file: Zephyr Biomedicals (A Division of Tulip Diagnostics Pvt. Ltd).

Important Note:

- The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.
- It is recommended to use the multiple channel pipettes to avoid time effect. A full plate of 96 wells may be used if automated pipetting is available.
- Duplication of standards & samples is not mandatory but may provide information on reproducibility & application errors.












LIMITATIONS OF THE ASSAY

- As with all diagnostic tests, a definite clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- The activity of the enzyme used is temperature-dependent and the RLU values may vary. The higher the room temperature (+18°C to +25°C) during substrate incubation, the greater will be the RLU values. Corresponding variations apply also to the incubation times. However, the Standards are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.
- Adaptation of this assay for use with automated sample processors and other liquid handling devices, in whole or in part, may yield differences in test results from those obtained using the manual procedure. It is the responsibility of each laboratory to validate that their automated procedure yields test results within acceptable limits.
- Insufficient washing (e.g., less than 5 wash cycles, too small wash buffer volumes, or shortened reaction times) can lead to incorrect results.

BIBLIOGRAPHY

- White, D.; Kramer, D.; Johnson, G.; Dick, F. And Hamilton, H. A.m. J. Clin. Path. 72: 346; 1986.
- Valberg, L. CMAJ. 122:1240; 1980.
- Forman, D. and Parker, S. Ann. Clin. Lab. Sci. 10: 345; 1980.
- Hazard, J.T.; Yokota, M.; Arosio, P. And Drysdale, J. Blood. 49:139; 1977.
- Smimes, M.A.; Addiego, Jr.J.E. and Dallman, P.R. Blood. 43:581; 1974.

SYMBOL KEYS

 Temperature Limitation	 Consult Instructions for use	 Date of Manufacture	 Batch Number / Lot Number
 Manufacturer	 In vitro Diagnostic Medical Device	 This side up	 Contains sufficient for <n> tests
 Use by	 Catalogue Number	 Do not reuse	

Manufactured by:
Zephyr Biomedicals

A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa- 403 722, INDIA.

Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.

An ISO 13485
Certified Company

0321/VER-01



Chemiluminescence Assay for the Quantitative Determination of Ferritin in Human Serum.
FOR IN VITRO DIAGNOSTIC USE ONLY
Store at 2°C to 8°C

INTENDED USE

ELECTRA™ Ferritin CLIA test is intended for the quantitative determination of Ferritin in human serum. For In Vitro Diagnostic Use only.

INTRODUCTION

One of the most prevalent disorders of man is the dietary deficiency of iron and the resulting anemia. Therefore, the assays of iron, total iron binding capacity and other assessments of iron compounds in the body are clinically significant.

Iron-storage compounds in the body include hemoglobin, hemosiderin, myoglobin and the cytochromes. In most tissues, ferritin is a major iron-storage protein. Human ferritin has a molecular weight of approximately 450,000 daltons, and consists of a protein shell around an iron core; each molecule of ferritin may contain as many as 4,000 iron atoms. Under normal conditions, this may represent 25% of the total iron found in the body.

High concentrations of ferritin are found in the cytoplasm of the reticuloendothelial system, the liver, spleen and bone marrow. Methods previously used to measure iron in such tissues are invasive, cause patient trauma and lack adequate sensitivity. The measurement of ferritin in serum is useful in determining changes in body iron storage, and is non-invasive with relatively little patient discomfort. Serum ferritin levels can be measured routinely and are particularly useful in the early detection of iron-deficiency anemia in apparently healthy people. Serum ferritin measurements are also clinically significant in the monitoring of the iron status of pregnant women, blood donors, and renal dialysis patients. High ferritin levels may indicate iron overload without apparent liver damage, as may be noted in the early stages of idiopathic hemochromatosis. Ferritin levels in serum have also been used to evaluate clinical conditions not related to iron storage, including inflammation, chronic liver disease, and malignancy.

PRINCIPLE

ELECTRA™ Ferritin Quantitative CLIA assay is for use on **ELECTRA** analyzers. **ELECTRA™ Ferritin CLIA** works on the principle of chemiluminescence wherein light is produced by a chemical reaction from a substance as it returns from an electronically excited state to the ground state. When catalysed by HRP, the oxidation of luminol by hydrogen peroxide produces an electronically excited form of 3-aminophthalate which on relaxation emits light with maximum intensity at $\lambda=425\text{nm}$.

The **ELECTRA™ Ferritin CLIA** system utilizes mouse monoclonal anti-ferritin antibody for solid phase (microwells) immobilization and another anti-ferritin antibody in the antibody-enzyme (horseradish peroxidase) conjugate solution. The test sample is allowed to react simultaneously with the two antibodies, resulting in the ferritin molecules being sandwiched between the solid phase and enzyme-linked antibodies. After incubation, the wells are washed with wash buffer to remove unbound anti-ferritin conjugate. A solution of chemiluminescent substrate is then added and Luminescence is measured in RLU. The intensity of the emitting light is proportional to the amount of enzyme present and is directly related to the amount of ferritin in the sample. By reference to a series of ferritin standards assayed in the same way, the concentration of ferritin in the unknown sample is quantified.

MATERIALS & COMPONENTS

Materials provided with the test kits:

- Coated Microwells: Microwells coated with anti-ferritin antibody.
- Enzyme Conjugate Reagent. Ready to use.
- Substrate A: Chemiluminescent substrate containing enhanced luminol solution.
- Substrate B: Chemiluminescent substrate containing stabilized peroxide solution.
- Ferritin Standard set of 6 standards labeled as A to F in liquid form. Ready to use. For standard Concentrations refer vial label.
- Wash Buffer Concentrate (20X).

Materials required but not provided:

- Precision pipettes: 10-100µl, 20-200µl, 100-1000µl
- Disposable pipette tips
- Distilled water
- Disposable Gloves
- ELECTRA™ Analyzer**



Chemiluminescence assay



Chemiluminescence assay



Chemiluminescence assay



Chemiluminescence assay

STORAGE AND STABILITY

1. **ELECTRA™ Ferritin** kit is stable at 2-8°C up to the expiry date printed on the label.
2. Coated microwells should be used within one month upon opening the pouch provided that once opened, the pouch must be resealed to protect from moisture. If the colour of the dessicant has changed from blue to white at the time of opening the pouch, another coated microwells pouch should be used.
3. Diluted Wash Buffer is stable upto one week when stored at 2-8°C.
4. Working Substrate (A+B) must be used immediately.

SPECIMEN COLLECTION

1. Collect blood specimen by venipuncture according to the standard procedure.
2. Only serum should be used.
3. Avoid grossly hemolytic, lipemic or turbid samples.
4. Preferably use fresh samples. However, specimens can be stored up to 48 hours at 2-8°C, for short duration.
5. For longer storage, specimens can be frozen at -20°C. Thawed samples must be mixed prior to testing.
6. Do not heat inactivate before use.
7. Specimen containing precipitate or particulate matter should be clarified by centrifugation prior to use.
8. Specimen should be free from particulate matter and microbial contamination.

PRECAUTIONS

1. Bring all reagents and specimen to room temperature before use.
2. Do not pipette any material by mouth.
3. Do not eat, drink or smoke in the area where testing is done.
4. Use protective clothing and wear gloves when handling samples.
5. Use absorbent sheet to cover the working area.
6. Immediately clean up any spills with sodium hypochlorite.
7. All specimens and standards should be considered potentially infectious and discarded appropriately.
8. Neutralize acid containing waste before adding hypochlorite.
9. Do not use kit after the expiry date.
10. Do not mix components of one kit with another.
11. Always use new tip for each specimen and reagent.
12. Do not allow liquid from one well to mix with other wells.
13. Do not let the strips dry in between the steps.


REAGENT PREPARATION

- All reagents should be brought to room temperature (18-25°C) and mixed by gently inverting or swirling prior to use. Do not induce foaming.
- Dilute Wash Buffer 20 times (for example add 5ml concentrated buffer to 95 ml distilled or deionized water). Mix well before use.
- Prepare a Working Substrate by Mixing Substrate A and Substrate B in equal volume (1:1 ratio) before addition to the micro-wells.

No. of Strips	1	2	3	4	5	6	7	8	9	10	11	12
Substrate-A µl	250	450	650	850	1050	1250	1450	1650	1850	2050	2250	2450
Substrate-B µl	250	450	650	850	1050	1250	1450	1650	1850	2050	2250	2450

TEST PROCEDURE

1. Secure the desired number of coated wells in the holder. Dispense **20µl** of Standards and Serums into the appropriate wells.
2. Dispense **100µl** of Enzyme Conjugate Reagent into each well. Incubate at room temperature(18-25°C) for **30 mins**.
3. Remove the incubation mixture by emptying the plate content into a waste container. Rinse and empty the microtiter plate 5 times with Wash Buffer (1X). Strike the microtiter plate sharply onto absorbent paper or paper towels to remove all residual water droplets.
4. Add **50µl** of working Substrate (A+B) in all the micro-wells. Keep away from direct light while adding the substrate.
5. Cover the **ELECTRA™** microplate and incubate for **10 minutes** at room temperature (**18-25°C**) in dark.
6. Read the **ELECTRA™** micro-plate exactly at **10 minutes** in **ELECTRA™ Analyzer**. If **ELECTRA™** micro-plate is not read between 10-15 minutes the test results should be considered as invalid.



CALCULATION OF RESULTS


Construct a Standard curve by plotting the mean RLU obtained from each reference standard against its concentration in ng/ml on the graph paper, with RLU values on the vertical or Y axis and concentrations on the horizontal or X axis. Use the RLU values for each specimen to determine the corresponding concentration of Ferritin in ng/ml from the standard curve. Any diluted specimens must be corrected by the appropriate dilution factor.

Example of Standard curve

Results of a typical Standard run with RLU's shown in the Y axis against Ferritin concentrations in the X axis.

Suggest: Use 4-Parameter Standard curve to calculate sample values.

Ferritin (ng/ml)	RLU's
A	71
B	19037
C	372175
D	871763
E	2854968
F	4549596



This Standard curve is for the purpose of illustration only, and should not be used to calculate samples. Each user should obtain his or her own Standard curve and data.

EXPECTED VALUES AND SENSITIVITY

Each laboratory must establish its own normal ranges based on patient population. The results provided below are from the literatures, which are based on a limited number of healthy adult blood specimens.

Adult Males	16-220 ng/ml
Adult Females	10-124 ng/ml
Newborn	22-220 ng/ml
Children (6 months – 15 years)	7-140 ng/ml

The minimum detectable concentration of Ferritin by this assay is estimated to be 5 ng/ml.