



Parameter	G-6PD
Intended Use	The Born Safe™ Neonatal G-6-PD Screening Assay is an enzymatic colorimetric assay for the quantitative determination of glucose-6-phosphate dehydrogenase (G6PD) activity in blood specimens dried on Whatman 903 filter paper. This test is intended for use as a screening method for red cell glucose-6-phosphate dehydrogenase deficiency in new-borns.
Principle	<p>The Neonatal G-6-PD Screening Assay uses dried blood spots specimens (cellulose paper) eluted in a buffer. After the elution step, the eluate containing G-6-PD is incubated with a reagent containing glucose-6-phosphate, which in the presence of NADP, catalyses the oxidation of glucose-6-phosphate to 6-phosphogluconate. The NADPH produced reacts with a colour reagent in which a tetrazolium salt gets reduced producing a distinct colour.</p> <p>This colour is measured at 550 nm and is directly proportional to the concentration of Glucose-6-phosphate dehydrogenase present in the sample. The results are calculated by evaluating the increase in OD per minute for unknowns against the calibrator with known G-6-PD activity. The reduced tetrazolium salt can be measured by a two-point measurement mode using two measurements one at time = 0 min, and the second one at time = 15 min.</p>
Kit Components	<p>Reagents</p> <ul style="list-style-type: none"> ● Calibrator and Controls blood spots: 1 + 1 set of dried blood spots cards of human whole blood spotted onto Whatman S&S 903 paper containing 1 calibrator and 2 controls. ● Elution Buffer: 1 X 8 ml of buffered solution. Ready to use. Preservative NaN_3 (<0.1%). ● Reagent: For 96T: 4 X 2 ml of Glucose-6-phosphate solution. Lyophilized. Reconstitute with 2 ml of distilled water. After reconstitution, the reagent can be stored at 2-8° C for one month. Preservative: NaN_3 (<0.1%). ● Colour Reagent: 1 X 8 ml. of tetrazolium salt. Ready to use. Preservative NaN_3 (<0.1%). ● Colour Booster: 1 X 1 ml of intermediate electron acceptor. Ready to use. Preservative NaN_3 (<0.1%).
Accessories	<ul style="list-style-type: none"> ● Round bottom microtiter plates (Elution Plates) ● Flat-bottom microtiter plates with superior optical quality (Assay Plates)
Pack size	96 Tests

Key Distinctions

🦶 Born Safe G6PD is Indian FDA approved, ISO certified. 🦶 Comprehensive External Clinical Evaluation, exclusively on neonatal population: >3200 samples. 🦶 Absolute Compliance with Lab QMS trend (e.g. CDC-PT, G6PD-PMF, Taiwan). 🦶 Born Safe G6PD Calibrators & assay kits are validated by CDC samples. 🦶 Excellent correlation to other commercial Colorimetric & Fluorometric assays and CDC target values. 🦶 Born Safe kits are available both in smaller (96T) and in larger pack sizes (192T) suitably fitting into the requirements of every NBS Lab. 🦶 Manufactured in India Plant, hence efficient management of supply & logistics across country. 🦶 **First of its kind:** Integral kit components e.g. Round bottom breakapart wells, Reaction wells are supplied along with other assay reagent components without any extra charge.