

# Labsafe™

**Disinfectant solution for laboratoryware**



**Safe Lab. Safe You.**

# Labsafe™

## Disinfectant Solution for Laboratory Ware

**LABSAFE** is a clear, pink red solution with a characteristic phenolic odour and a broad-spectrum microbicidal activity. **LABSAFE** is recommended for decontamination of laboratory ware prior to disposal. Its use ensures complete safety to laboratory personnel.



### Applications:

#### MEDICAL/RESEARCH LABS

Heat sensitive laboratory ware  
Laboratory ware prior to disposal  
Pre-treatment of infectious waste  
TB reference labs

#### INDUSTRY

R & D and QA/QC labs of  
pharmaceutical, food and dairy industry

### Benefits :

- Broad spectrum biocidal action
- Pronounced tuberculocidal action
- Sporicidal in just 16 hours
- Effective decontaminant action
- Long shelf-life
- Good material compatibility


*“Phenol has mycobactericidal activity...”<sup>1, 2, 3</sup>*

*“Phenols most effective against microorganisms are moderately affected by organic matter”<sup>4</sup>*

#### REFERENCES:

1. Rubin, J. 1983. Agents for disinfection and control of tuberculosis. In disinfection, sterilization and preservation, 3<sup>rd</sup> edn(ed. Block, S S), Philadelphia: Lea@ febiger. 414-421
2. Gardner, J E and M M Peel. 1986. Introduction to sterilization and disinfection..Edinburgh:Churchill Livingstone.
3. Favero, M.S.1991. Practical application of liquid sterilants in health care facilities. In Sterilisation of medical products.(eds. Morrissey,R.F.and propopenko, Y I)volume V.
4. Key characteristics of the principal hard surface disinfectants. National Trade Publications, Inc. <http://www.corpchem.com/PDF/Advice/KeyCharactDisinf.pdf>

Distributed by:



*Applying Science In Disinfection*

Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403202, INDIA  
Tel.: +91 832 2458546-50. Fax: +91 832 2458544 E-mail: [sales@tulipgroup.com](mailto:sales@tulipgroup.com). Website: [www.tulipgroup.com](http://www.tulipgroup.com)



For Professional use only.  
Refer label before use.

**Safe Lab. Safe You.**

LSE/5000/06-09/----/V1