

**Graviscreeen<sup>®</sup>**  
**Rapid one step test dipstick for hCG**

**Intended Use:** \_\_\_\_\_

Graviscreeen<sup>®</sup>, a rapid one step test dipstick for the detection of hCG in urine during pregnancy.

**Graviscreeen<sup>®</sup> Dipstick kit components:**

<b>Dipstick</b>	Membrane assembly pre-dispensed with monoclonal Anti-hCG antisera - colloidal gold conjugate, rabbit IgG colloidal gold conjugate monoclonal Anti-hCG antiserum at 'Test region' and anti-rabbit anti serum at the 'Control region'
<b>Pipette</b>	Disposable plastic sample dropper
<b>Other Accessories</b>	Desiccant pouch and Package Insert

**Performance Characteristics:**

<b>Sensitivity:</b> ~ 20 mIU/ml of hCG	<b>Specificity:</b> 100% specific to $\beta$ hCG
--	--

FEATURES	BENEFITS
<b>Sensitivity adjusted to detect &gt; 20 mIU/ml of hCG</b>	Optimum sensitivity for safe and confident reporting of results Facilitates early pregnancy detection
<b>Specific for beta sub unit of hCG molecule</b>	No cross-reactivity with other hormones
<b>Storage: 4-30°C</b>	Flexibility of remote distribution and on the spot testing
<b>Non invasive procedure</b>	Routine pregnancy testing using urine specimens High patient compliance
<b>Composite kit</b>	No additional material required for testing

Storage / Stability	Temperature	Duration
<b>Unopened kit</b>	4-30°C	24 Months
<b>Opened kit</b>	4-30°C	< 1 Hour

Available Pack Sizes
25 Tests