

OCCULT BLOOD DETECTION



OCCULT BLOOD

Turnips, Cabbage
Watermelons, O
Red Meat, Cabb
Watermelons,
Broccoli, Cabb
Watermelons,
Turnips, Cabbage
Cabbage, Occult

wer, Beets, Radishes, Beets,
anges, Cauliflower, Tangelos,
Beets, Radishes, Broccoli,
anges, Vitamin C, Carrot,
ets, Radishes, Mushrooms,
anges, Beets, Cucumber,
eets, Radishes, Mushrooms,
anges, Vitamin C, Cabbage

Now, so easy so specific...

CANCHECK[®] - FOBT

Rapid test for detection of Faecal Occult Blood in human faeces

Presence of very small amount of blood (Occult blood) in stool is an early sign of colorectal bleeding or intestinal disorder such as colitis, polyps in colon or rectum and colorectal carcinomas (CRC). Conventional Gum - guaiac test or Benzidine based test for occult blood detection has limitations of false positive results due to dietary interference (interference due to food intake like cabbage, beet, cauliflower etc.) as well as false negative result in early stage of disease due to lack of sensitivity. **CANCHECK - FOBT** has overcome these limitations and provides reliable result through a simple user friendly test procedure.

High sensitivity, can detect colorectal bleeding in early stage - Detects occult blood as low as 200µg /Ltr. faeces.


Excellent specificity - Directly detect human hemoglobin to aid the diagnosis of colorectal bleeding. No cross reactivity of peroxidase reaction.

No food restriction for patients prior to the testing - No dietary interference in results.


Easy, user friendly test procedure and fast result - Result in 5 minutes, can be used as point of care testing.

Convenient pack size & storage - 10 tests pack, 4-30° C storage.


TEST PROCEDURE




INTERPRETATION OF RESULTS




NEGATIVE



POSITIVE



INVALID - Repeat the test



Early & specific detection of occult blood... **CANCHECK[®] - FOBT**

ISO 13485
Certified Company



TULIP DIAGNOSTICS (P) LTD.

For further information contact :

Gitanjali, Tulip Block, Dr. Antonio do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403 202, INDIA.

Tel.: +91 832 2458546-50 E-mail : sales@tulipgroup.com Website : www.tulipgroup.com