



	ANTI-K
Application	Monoclonal Anti-K (KELL) antibodies for slide and tube method
Titre	• ~ 1:32
Class of Antibody	IgM, Monoclonals
Source	Human cell line
Specificity	• ANTI-K 100% to K antigen
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation		Pack	Cat.No.
Rh phenotyping reagent(ANTI-K)	ANTI-K	2 ml	101750002
		5 ml	101750005
		10 ml	101750010

