



	<h2>Anti-D IgG</h2>
Application	Monoclonal Anti-D IgG reagent in a standardized concentration for weak D testing on Matrix™ Coombs Anti-IgG card
Reagent	Matrix™ Anti-D IgG contains: <ul style="list-style-type: none"> • Monoclonal Anti-D IgG reagent in a standardized concentration
Standardisation	Follows FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
Matrix™ Anti-D IgG	5 ml	102790005

