



APTT Reagent

Application	For the Activated Partial Thromboplastin Time determination in plasma, monitoring of heparin therapy and assessment of Intrinsic pathway
Reagent	Ready to use Liquid stable, Soyabean cephaloplastin reagent with Ellagic acid as an activator
Source	Soyabean phospholipid
Normal Range	23 - 33 seconds
Storage / Stability	2-8°C, 18 months 18-25°C, 1 week 37°C, 2 days

Presentation	Pack	Cat.No.
CELIN-SE	12 x 3 ml	10635123