



Chemiluminescence assay for the detection of Malaria specific antigen (pLDH) in human whole blood

Application	For detection of malaria specific antigen (pLDH) in human whole blood samples.
Principle	Chemiluminescence (CLIA)
Reagents	<ul style="list-style-type: none"> • Chemiluminescence compatible microwells coated with agglutinating sera for Pan malaria specific pLDH. • Controls for assay validation • Conjugate for antigen & antibody detection • Antibody Reagent: agglutinating sera for Pan malaria specific pLDH • Chemiluminescent substrate • Wash buffer
Storage / Stability	2-8°C, 12 months

Presentation	Pack	Cat.No.
Electra MAL Ag CLIA kit contains: • Microwell strips	96 Tests	407050096
	192 Tests	407050192

