



	Chemiluminescence assay for the detection of antibodies to <i>Treponema Pallidum</i> in human serum or plasma
Application	For detection of total antibodies (i.e IgG, IgA & IgM etc) to <i>Treponema Pallidum</i> in human serum or plasma.
Principle	Chemiluminescence (CLIA)
Reagents	<ul style="list-style-type: none"> • Chemiluminescence compatible microwells coated with recombinant 47Kd & 17Kd antigens • Controls for assay validation • Conjugate for antibody detection • Chemiluminescent substrate • Wash buffer
Storage / Stability	2-8°C, 12 months

Presentation	Pack	Cat.No.
Electra Tp. Ab. CLIA kit contains: • Microwell strips	96 Tests	407040096
	192 Tests	407040192

