



<p>FIAcheck™ FSH</p>	<p>Fluorescence Immunoassay for quantitative determination of FSH (Follicle stimulating hormone) in human serum</p>
<p>Application</p>	<p>FIAcheck FSH Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of FSH (Follicle stimulating hormone) in human serum.</p>
<p>Principle</p>	<p>Rapid, lateral flow, quantitative immuno-fluorescence assay</p>
<p>Measuring range</p>	<p>1-200mIU/mL</p>
<p>Reagents</p>	<p>Membrane assembly coated with</p> <ul style="list-style-type: none"> •Immobilized FSH mouse monoclonal antibody1 on test line •Goat anti chicken IgY in control line •Mixture of FSH mouse monoclonal antibody2 and chicken IgY labeled with fluorescent microspheres on the binding pad <p>Sample diluent</p>
<p>Test duration</p>	<p>15 mins</p>
<p>Storage / Stability</p>	<p>4-30°C, 12 months</p>

Presentation	Pack	Cat.No.
<p>FIAcheck FSH kit contains an individually pouched device and calibration QR code card</p>	<p>25 Tests</p>	<p>535010025</p>