



BabySafeTM

Newborn Screening

Key Features:

- Light Reflection based Rapid NBS test Card Reader
- Lot Specific Calibration through QR
- Quick, Simple & user-friendly software
- Portable and Compact: Easy to carry and use anywhere
- Faster Turnaround Time



Technical Specifications

| | |
|------------------------------|--|
| Principle | Reflection of light |
| Display | 7-inch display |
| Language | English |
| Port | USB: 2 Nos LAN: 2 Nos RS232: 1 Nos end |
| Clock | Inbuilt RTC |
| Individual Test Speed | 30 seconds |
| Result | Quantitative / Qualitative Test |
| Printer | Inbuilt Thermal printer |
| Paper Roll | 56 mm paper roll |
| Power Supply | DC 24V, 2.5Amp External Power adaptor |
| Dimension (L×W×H) mm | 265×250×125 mm |
| Weight | 2.5Kg |

The **BabysafeTM** analyzer is providing in a single carton with accessories described below:

| Components | Qty. | Components | Qty. |
|---------------------------------|------|--|-------|
| Babysafe TM Analyzer | 1 | User manual | 1 |
| Power Adapter | 1 | Pipette (2 μ l, 10 μ l & 50 μ l) | 1 set |
| Thermal Paper Roll | 1 | | |



TULIP DIAGNOSTICS (P) LTD.