

Probe Cleaner VET

(Probe Cleaner For Hematology Analyzer)

(For veterinary invitro diagnostic use only)

INTENDED USE

Probe Cleaner VET is used on Hematology Analyzer for cleaning & washing tubing system, removing blood residual & other particles to ensure blood cell counting. Probe Cleaner removes the cell residues, blood proteins in the instrument to prevent protein precipitation in the sampling needle, pipe, counting chamber and prolongs the life of the machine, and helps in the smooth functioning of the system.

PRESENTATION

REF	1126270050
Pack Size	50 ml
Probe Cleaner VET	50 ml

COMPOSITION

Sodium Hypochlorite: > 100 mMol.

STORAGE & STABILITY

This reagent should be stored in a clean dry area at 23-29°C. Keep the product away from direct sunlight.

The product is stable till the expiry mentioned on the label.

Once opened, the reagent loses active chlorine hence, close the cap tightly immediately after use.

PROCEDURE

Refer to the instrument manual & the program given there in for the relevant instrument.

CAUTIONS

1. Do not pipette by mouth or inhale the reagent.
2. If the reagent comes in contact with eye or skin, wash off immediately with plenty of water & seek medical attention if necessary.
3. Please do not use this reagent if there are any precipitates, particles or floccules in it.
4. Please dispose the waste & contaminated packing material according to local laws & regulations.

SYMBOL KEYS

 Store at 23-29°C	 Manufacturer	 <i>In vitro</i> Diagnostic Medical Device	 Date of Manufacture	 This way up
 Use by (Last day of stated month)	 Consult Instructions for use	 Batch Number	 Catalogue Number	

CCLV10625/VER-01



Manufactured by:

Coral Clinical Systems

A Division of Tulip Diagnostics (P) Ltd.

BUILDING E, PLOT NO. M-46/47, PHASE III B, VERNA INDUSTRIAL ESTATE, VERNA, GOA-403 722, INDIA.

REGD. OFFICE : GITANJALI, TULIP BLOCK, DR. ANTONIO DO REGO BAGH,

ALTO SANTACRUZ, BAMBOLIM COMPLEX P.O., GOA-403 202, INDIA.