



QUADRAPEDTM
Mastindex CMT Reagent

Mastindex California Mastitis Test Reagent

Application	Somatic Cell Count in milk samples from individual animals can be performed using Mastindex CMT Reagent
Reagent	Contains ready to use anionic surfactant with bromocresol purple
Principle	When milk and Mastindex CMT Reagent is mixed in equal amount, the reagent dissolves and forms a stringy mass in presence of clinically significant somatic cells and homogenous mixture in the absence of clinically significant somatic cells
Test duration	One-minute assay
Storage/ Stability	Room temperature / 24 months

Presentation	Pack Size	Cat. No.
QUADRAPEDTM Mastindex CMT Reagent	500 ml	204170010500