



 ALBUMIN KIT	ALBUMIN KIT (BCG Method)
Application	For the determination of Albumin in serum or plasma
Reagent	Ready to use
Storage/ Stability	Reagent :23 - 29°C/ 24 months Standard : 2 - 8°C/ 24 months
Specific Performance Characteristics	
LOB	0.024 g/dl at 630 nm
LOD	0.1 g/dl at 630 nm
LOQ	0.25 g/dl at 630 nm
Lower Limit	0.1 g/dl
Higher Limit	7 g/dl

Presentation	Pack Size	Cat. No.
QUADRAPED ALBUMIN KIT	25 ml	1126010025