



<p style="text-align: center;">  ALKALINE PHOSPHATASE KIT (DEA) </p>	<p style="text-align: center;"> ALKALINE PHOSPHATASE KIT (DEA) (pNPP Kinetic Method) </p>
<p>Application</p>	<p>For the determination of Alkaline Phosphatase Activity in serum.</p>
<p>Reagent</p>	<p>Substrate Tablets and Buffer Reagent Working Reagent preparation: Dissolve 1 Substrate Tablet in 3.2 ml of Buffer Reagent</p>
<p>Specific Performance Characteristics</p>	
<p>LOD</p>	<p>3 U/L</p>
<p>LOQ</p>	<p>28 U/L</p>
<p>Lower Limit</p>	<p>3 U/L</p>
<p>Higher Limit</p>	<p>700 U/L</p>
<p>Storage/ Stability</p>	<p>2-8°C / 18 months</p>

Presentation	Pack Size	Cat. No.
QUADRAPED ALKALINE PHOSPHATASE KIT (DEA)	10 x 3 ml	1126020103