



<p style="text-align: center;">QUADRAPEDTM BILIRUBIN KIT</p>	<p style="text-align: center;">BILIRUBIN KIT (Mod. Jendrassik and Grof's Method)</p>
<p>Application</p>	<p>For the determination of Direct and Total Bilirubin in serum.</p>
<p>Reagent</p>	<p>Ready to use</p>
<p>Specific Performance Characteristics</p>	
<p>LOD</p>	<p>0.03 mg/dl</p>
<p>LOQ</p>	<p>0.05 mg/dl</p>
<p>Lower Limit</p>	<p>0.03 mg/dl</p>
<p>Higher Limit</p>	<p>20 mg/dl</p>
<p>Storage/ Stability</p>	<p>Reagent :23-29°C/ 24 months Standard : 2-8°C/ 24 months</p>

Presentation	Pack Size	Cat. No.
QUADRAPED BILIRUBIN KIT	35 Tests	1126040035