



QUADRAPED™
BRUCEL-RB

**Slide Screening Test For
Brucella Antibodies**

Application	Rapid slide agglutination test for screening of antibodies to Brucella in animal serum including cattle, sheep, goats, and swine
Specificity	Specific reactivity towards antibodies to Brucella antigens
Reagent	Smooth, killed buffered suspensions of Brucella abortus strain 99, coloured with Rose Bengal
Positive control	Reactive with the reagent
Negative control	Non-Reactive with the reagent
Test duration	One-minute assay
Storage/ Stability	2-8°C/ 24 months

Presentation	Pack Size	Cat. No.
QUADRAPED™ BRUCELLA-RB	100 Tests	105960100