



**QUADRAPED<sup>TM</sup>**  
**CALCIUM KIT**

**CALCIUM KIT**  
**(Arsenazo III Method)**

<b>Application</b>	For the determination of Calcium in serum or plasma.
<b>Reagent</b>	Ready to use
<b>Specific Performance Characteristics</b>	
LOD	0.05 mg/dl at 630nm
LOQ	0.625 mg/dl at 630nm
Lower Limit	0.274 mg/dl
Higher Limit	15 mg/dl
<b>Storage/ Stability</b>	2 - 8°C/ 24 months

<b>Presentation</b>	<b>Pack Size</b>	<b>Cat. No.</b>
<b>QUADRAPED CALCIUM KIT</b>	25 tests	1126050025