



QUADRAPED™
CHOLESTEROL SR KIT

**CHOLESTEROL SR KIT
(CHOD / PAP Method)**

| | |
|---|--|
| Application | For the determination of Cholesterol in serum or plasma. |
| Reagent | Ready to use |
| Specific Performance Characteristics | |
| LOD | 1 mg/dl |
| LOQ | 12.5 mg/dl |
| LOB | 4.8 mg/dl |
| Lower Limit | 6.25 mg/dl |
| Higher Limit | 750 mg/dl |
| Storage/ Stability | 2 - 8°C/ 18 months |

| | | |
|------------------------------|------------------|-----------------|
| Presentation | Pack Size | Cat. No. |
| QUADRAPED CHOLESTEROL SR KIT | 25 ml | 1126060025 |