



 CREATININE KIT	CREATININE KIT (Mod. Jaffe's Kinetic Method)
Application	For the determination of Creatinine in serum
Reagent	Ready to use
Specific Performance Characteristics	
LOD	0.1 mg/dl
LOQ	0.06 mg/dl
Lower Limit	0.1 mg/dl
Higher Limit	20 mg/dl
Storage/ Stability	23 - 29°C/ 24 months

Presentation	Pack Size	Cat. No.
QUADRAPED CREATININE KIT	2 x 35 ml	1126080235