



<p style="text-align: center;">   <b>ELYTE 3 KIT</b> </p>	<p style="text-align: center;"><b>ELYTE 3 KIT (Colorimetric Method)</b></p>
<p><b>Application</b></p>	<p>For the determination of sodium, potassium &amp; chloride in serum.</p>
<p><b>Reagent</b></p>	<p>Ready to use</p>
<p><b>Specific Performance Characteristics</b></p>	
<p><b>Sodium</b></p>	
<p>LOD</p>	<p>1.0 mmol/l</p>
<p>LOQ</p>	<p>5.0 mmol/l</p>
<p>Lower Limit</p>	<p>1.0 mmol/l</p>
<p>Higher Limit</p>	<p>200 mmol/l</p>
<p><b>Potassium</b></p>	
<p>LOB</p>	<p>0.07 mmol/l</p>
<p>LOQ</p>	<p>0.31 mmol/l</p>
<p>LOD</p>	<p>0.15 mmol/l</p>
<p>Lower Limit</p>	<p>0.15 mmol/l</p>
<p>Higher Limit</p>	<p>8 mmol/l</p>
<p><b>Chloride</b></p>	
<p>LOB</p>	<p>0.29 mmol/l</p>
<p>LOQ</p>	<p>5.0 mmol/l</p>
<p>LOD</p>	<p>1.0 mmol/l</p>
<p>Lower Limit</p>	<p>1.0 mmol/l</p>
<p>Higher Limit</p>	<p>140 mmol/l</p>

<b>Storage/ Stability</b>	23-29°C/ 24 months
---------------------------	--------------------

<b>Presentation</b>	<b>Pack Size</b>	<b>Cat. No.</b>
<b>QUADRAPED ELYTE 3 KIT</b>	15 Tests	1126090015