



	<b>GGT KIT (Carboxy Substrate Method)</b>
<b>Application</b>	For the determination of g Glutamyl Transferase Activity in serum
<b>Reagent</b>	Substrate Tablets and Buffer Reagent Working Reagent preparation: Dissolve 1 Substrate Tablet in 2.2 ml of Buffer Reagent
<b>Specific Performance Characteristics</b>	
LOD	2 U/L
LOQ	4 U/L
Lower Limit	2 U/L
Higher Limit	700 U/L
<b>Storage/ Stability</b>	2-8°C/ 18 months Working reagent 15 days at 2-8°C

Presentation	Pack Size	Cat. No.
QUADRAPED GGT KIT	10 x 2 ml	1126110102