



 QUADRAPEDTM GLUCOSE KIT	GLUCOSE KIT (GOD / POD Method)
Application	For the determination of Glucose in serum or plasma.
Reagent	Ready to use
Specific Performance Characteristics	
LOB	1.42 mg/dl at 505 nm
LOD	1.90 mg/dl at 505 nm
LOQ	3.13 mg/dl at 505 nm
Lower Limit	1.90 mg/dl
Higher Limit	500 mg/dl
Storage/ Stability	2 - 8°C/ 24 months

Presentation	Pack Size	Cat. No.
QUADRAPED GLUCOSE KIT	50 ml	1126100050