



<p style="text-align: center;">QUADRAPEDTM HDL-D CHOLESTEROL KIT</p>	<p style="text-align: center;">HDL-D CHOLESTEROL KIT (Direct Enzymatic Method)</p>
<p>Application</p>	<p>For the determination of HDL Cholesterol in serum/plasma.</p>
<p>Reagent</p>	<p>Ready to use</p>
<p>Specific Performance Characteristics</p>	
<p>LOD</p>	<p>1.0 mg/dl</p>
<p>LOQ</p>	<p>2.0 mg/dl</p>
<p>Lower Limit</p>	<p>1.0 mg/dl</p>
<p>Higher Limit</p>	<p>150 mg/dl</p>
<p>Storage/ Stability</p>	<p>2 - 8°C/ 18 months</p>

Presentation	Pack Size	Cat. No.
QUADRAPED HDL-D CHOLESTEROL KIT	40 ml	1126120040