



 LDH (P-L) KIT	LDH (P-L) KIT (DGKC METHOD)
Application	For the determination of Lactate Dehydrogenase Activity in serum.
Reagent	Working reagent is prepared by mixing Starter (L2) reagent and Buffer (L1) reagent
Specific Performance Characteristics	
LOD	8 U/L
LOQ	10 U/L
Lower Limit	8 U/L
Higher Limit	2000 U/L
Storage/ Stability	2-8°C / 18 months Working reagent 2-8°C/ 3 weeks

Presentation	Pack Size	Cat. No.
QUADRAPED LDH (P-L) KIT	25 ml	1126130025