



<p style="text-align: center;">QUADRAPEDTM LDL-D CHOLESTEROL KIT</p>	<p style="text-align: center;">LDL-D CHOLESTEROL KIT (Direct Enzymatic Method)</p>
<p>Application</p>	<p>For the determination of LDL Cholesterol in serum / plasma.</p>
<p>Reagent</p>	<p>Ready to use</p>
<p>Specific Performance Characteristics</p>	
<p>LOD</p>	<p>2 mg/dl</p>
<p>LOQ</p>	<p>6 mg/dl</p>
<p>Lower Limit</p>	<p>2 mg/dl</p>
<p>Higher Limit</p>	<p>1000 mg/dl</p>
<p>Storage/ Stability</p>	<p>2 - 8°C/ 18 months</p>

<p>Presentation</p>	<p>Pack Size</p>	<p>Cat. No.</p>
<p>QUADRAPED LDL-D CHOLESTEROL KIT</p>	<p>40 ml</p>	<p>1126140040</p>