



 QUADRAPEDTM LIPASE KIT	LIPASE KIT (Colorimetric Method)
Application	For determination of Lipase Activity in serum / plasma.
Reagent	Ready to use
Specific Performance Characteristics	
LOD	2.5 U/L
LOQ	5 U/L
Lower Limit	2.5 U/L
Higher Limit	300 U/L
Storage/ Stability	2-8°C / 18 months

Presentation	Pack Size	Cat. No.
QUADRAPED LIPASE KIT	10 ml	1126150010