



 <b>QUADRAPED™</b> <b>MAGNESIUM KIT</b>	<b>MAGNESIUM KIT</b> <b>(CALMAGITE METHOD)</b>
<b>Application</b>	For the determination of Magnesium in serum.
<b>Reagent</b>	Ready to use
<b>Specific Performance Characteristics</b>	
LOD	0.05 mEq/L
LOQ	0.55 mEq/L
Lower Limit	0.05 mEq/L
Higher Limit	6.0 mEq/L
<b>Storage/ Stability</b>	2 - 8°C/ 18 months

<b>Presentation</b>	<b>Pack Size</b>	<b>Cat. No.</b>
QUADRAPED MAGNESIUM KIT	25 ml	1126160025