



	PHOSPHORUS KIT (Molybdate U.V. Method)
Application	For the determination of Inorganic Phosphorus in serum or plasma.
Reagent	Working reagent is prepared by mixing Molybdate (L2) reagent and Acid (L1) reagent
Specific Performance Characteristics	
LOB	0.1 mg/dl
LOD	0.1 mg/dl
LOQ	0.62 mg/dl
Lower Limit	0.1 mg/dl
Higher Limit	20 mg/dl
Storage/ Stability	23-29°C/ 24 months Working reagent 2-8°C/ 6 months

Presentation	Pack Size	Cat. No.
QUADRAPED PHOSPHORUS KIT	25 ml	1126170025