



 <b>SGOT KIT</b>	<b>SGOT (ASAT) KIT (Mod. IFCC Method)</b>
<b>Application</b>	For determination of SGOT (ASAT) Activity in serum.
<b>Reagent</b>	Working reagent is prepared by mixing Starter (L2) reagent and Enzyme (L1) reagent
<b>Specific Performance Characteristics</b>	
LOD	8.6 U/L
LOQ	17.1 U/L
Lower Limit	8.6 U/L
Higher Limit	500 U/L
<b>Storage/ Stability</b>	2-8°C / 18 months Working reagent 2-8°C/ 3 weeks

<b>Presentation</b>	<b>Pack Size</b>	<b>Cat. No.</b>
QUADRAPED SGOT KIT	25 ml	1126180025