



 SGPT KIT	SGPT (ALAT) KIT (Mod. IFCC Method)
Application	For determination of SGOT (ASAT) Activity in serum.
Reagent	Working reagent is prepared by mixing Starter (L2) reagent and Enzyme (L1) reagent
Specific Performance Characteristics	
LOD	3.8 U/L
LOQ	7.5 U/L
Lower Limit	3.8 U/L
Higher Limit	500 U/L
Storage/ Stability	2-8°C / 18 months Working reagent 2-8°C/ 3 weeks

Presentation	Pack Size	Cat. No.
QUADRAPED SGPT KIT	25 ml	1126190025