



 TOTAL PROTEIN KIT	TOTAL PROTEIN KIT (Biuret Method)
Application	For determination of Total Proteins in serum and plasma.
Reagent	Ready to use
Specific Performance Characteristics	
LOB	0.079 g/dl at 546 nm
LOD	0.1 g/dl at 546 nm
LOQ	1.0 g/dl at 546 nm
Lower Limit	0.1 g/dl
Higher Limit	15 g/dl
Storage/ Stability	2-8°C / 24 months

Presentation	Pack Size	Cat. No.
QUADRAPED TOTAL PROTEIN KIT	25 ml	1126200025