



 TRIGLYCERIDES SR KIT	TRIGLYCERIDES SR KIT (GPO / PAP Method)
Application	For determination of Triglycerides in serum or plasma.
Reagent	Ready to use
Specific Performance Characteristics	
LOD	1.0 mg/dl
LOQ	5.0 mg/dl
Lower Limit	1.0 mg/dl
Higher Limit	1000 mg/dl
Storage/ Stability	2-8°C / 18 months

Presentation	Pack Size	Cat. No.
QUADRAPED TRIGLYCERIDES SR KIT	25 ml	1126210025