



| | |
|---|--|
|  URIC ACID KIT | URIC ACID KIT (Uricase / PAP Method) |
| Application | For determination of Uric Acid in serum or plasma. |
| Reagent | Ready to use |
| Specific Performance Characteristics | |
| LOD | 0.15 mg/dl |
| LOQ | 1.0 mg/dl |
| Lower Limit | 0.15 mg/dl |
| Higher Limit | 20 mg/dl |
| Storage/ Stability | 2-8°C / 24 months |

| | | |
|--------------------------------|------------------|-----------------|
| Presentation | Pack Size | Cat. No. |
| QUADRAPED URIC ACID KIT | 25 ml | 1126230025 |