



 <b>ZINC KIT</b>	<b>ZINC KIT (Colorimetric Method)</b>
<b>Application</b>	For determination of Zinc in serum.
<b>Reagent</b>	Ready to use
<b>Specific Performance Characteristics</b>	
LOD	5.0 µg/dl
LOQ	10.0 µg/dl
Lower Limit	5.0 µg/dl
Higher Limit	700.0 µg/dl
<b>Storage/ Stability</b>	2-8°C / 18 months Working reagent 2-8°C / 3 weeks

<b>Presentation</b>	<b>Pack Size</b>	<b>Cat. No.</b>
QUADRAPED ZINC KIT	25 ml	1126240025