



CombiVIEW™ Pv/Pf	Rapid test for malaria Pv/Pf (Pf HRP2 + Pf pLDH and Pv pLDH)
Application	For the detection of <i>P. falciparum</i> specific HRP 2 + <i>p. falciparum</i> specific pLDH and <i>P. vivax</i> specific pLDH in whole blood
Principle	Two site sandwich immunoassay, lateral flow device
Sensitivity	100% correlation with microscopy
Specificity	100% specific to <i>P. falciparum</i> HRP 2, <i>p. falciparum</i> pLDH and <i>P. vivax</i> specific pLDH
Reagents	Membrane assembly pre-dispensed with • Agglutinating Sera for Pf HRP-2 + Pf pLDH - colloidal gold conjugate • Agglutinating Sera for <i>P. vivax</i> specific pLDH colloidal gold • Agglutinating Sera for Pf HRP-2+ Pf pLDH • Agglutinating Sera for Malaria specific pLDH • Agglutinating Sera for Rabbit globulin at the respective regions. • Clearing Buffer
Test duration	Twenty minutes assay
Storage / Stability	4-40°C, 24 months

Presentation	Pack	Cat.No.
Combiview kit contains individually pouched device with clearing buffer and 5 µl sample applicator	10 T 25 T 25 X 1T	503230010 503230025 503230025 (1T)