



	Rapid test for <i>P. falciparum</i> malaria (Pf HRP2 and Pf pLDH)
Application	For the differential detection of <i>P. falciparum</i> specific HRP 2 in whole blood
Principle	Two site sandwich immunoassay, lateral flow device
Sensitivity	> 99% correlation with microscopy
Specificity	100% specific to <i>P.falciparum</i> HRP 2 and <i>p. falciparum</i> pLDH
Reagents	Membrane assembly predisposed with - • Agglutinating Sera for HRP-2 Pf pLDH • Colloidal gold conjugate, rabbit IgG colloidal gold conjugate – • Agglutinating Sera for HRP-2 at 'Test' region H • Agglutinating Sera for Pf pLDH at test region L • Agglutinating Sera for rabbit globulin at 'Control' region • Clearing Buffer
Test duration	Twenty minutes assay
Storage / Stability	4-40°C, 24 months

Presentation	Pack	Cat.No.
DuoCheck Pf kit contains individually pouched device with clearing buffer and 5 µl sample applicator	10 T	302110010
	25 T	302110025
	50 T	302110050
	100 T	302110100
	25 X 1T	302110025 (1T)