



	<b>Rapid Competitive Immunochromatographic Test for detection of amphetamine in urine</b>
<b>Application</b>	For the detection of presence of amphetamine in human urine specimens qualitatively with concentrations as low as 1000 ng/ml
<b>Principle</b>	Self performing, two site sandwich immunoassay, lateral flow device
<b>Reagents</b>	<p>A. Each INSIGHT AMP kit contains individual pouches each containing a</p> <ol style="list-style-type: none"> <li>1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti amphetamine antibody monoclonal and mouse IgG, amphetamine conjugated to BSA and goat anti mouse IgG at the respective regions.</li> <li>2. Sample dropper.</li> <li>3. Desiccant pouch.</li> </ol> <p>B. Package insert.</p>

Presentation	Pack	Cat.No.
<b>Insight™ AMP kit contains individually pouched device</b>	10 Tests	10803010

