



	Rapid One Step Test (device) for detection of HCG
Application	For the detection of human chorionic gonadotrophin (HCG) hormone in urine
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Sensitivity	~ 25 mIU/ml
Reagents	Each kit contains: A. Individual pouches, each containing a <ol style="list-style-type: none"> 1. Membrane assembly impregnated with Anti HCG antiserum- colloidal gold conjugate, rabbit IgG colloidal gold conjugate, anti-HCG antiserum and anti rabbit antiserum coated at the respective regions. 2. Disposable Plastic Dropper. 3. Desiccant pouch. B. Package insert.
Test duration	5 mins.
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
Tina™ HCG kit contains individually pouched device	50 Tests	707030050

