

**Paratrol™****Control set for immunological assay of *P. falciparum*, *P. vivax* Malaria antigen detection**

Intended Use: _____

Controls for on site validation / quality assurance of immunological assays for malaria antigen detection.

Paratrol™ Kit components:

Reagent	Set of lyophilized preparation of human whole blood sources with <i>P.falciparum</i> , <i>P. vivax</i> and malaria negative whole blood
Other Accessories	Package Insert

FEATURES	BENEFITS
Precisely prepared reagent-with known parasitic load	Appropriately reconstituted controls, ensures quality and standardized results of malaria antigen detection immuno assay
Freeze dried reagent	Aids in minimizing microbial contamination
Convenient pack size	Ideal for all users

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	12 Months
Reconstituted vial	2-8°C	2 Days

Available Pack Sizes
2 x 3 x 0.2 ml